

# Medical Device Alert

## Immediate action

Ref: MDA/2009/060 Issued: 22 September 2009 at 14:00

### Device

Dialysis bloodlines.  
AV-SET Dialog® GB  
Product code 7210697.  
Manufactured by BBraun Avitum AG.



### Problem

Blood leakage at the arterial chamber (PBE chamber) junction.

### Action by

All staff working in renal units and other areas where haemodialysis machines are used.

### CAS deadlines

Action underway: 06 October 2009  
Action complete: 13 October 2009

### Action

- Check lot numbers to identify whether you have any affected devices.
- Quarantine and do not use affected devices.
- Contact BBraun for replacements.
- If replacements are unavailable from BBraun, source sets from alternative suppliers.

### Contact

#### Manufacturer

Ms Catherine Clulow  
BBraun Medical Ltd  
Tel: 0114 225 9155  
Fax: 0114 225 9111

E-mail: [catherine.clulow@bbraun.com](mailto:catherine.clulow@bbraun.com)

## Link to full Medical Device Alert

<http://www.mhra.gov.uk/Publications/Safetywarnings/MedicalDeviceAlerts/CON057334>