

# Medical Device Alert

## Immediate action

Ref: MDA/2009/075 Issued: 18 November 2009 at 11:00

### Device

Needle-free intravascular connectors:

- BD Q-Syte stand-alone Luer access split septum connector
- BD Q-Syte with 15cm micro bore extension set.



Manufactured by Becton Dickinson.

NHS Supply Chain code FSC045 and FSC055 respectively.

Problem	Action
<p>The manufacturer has received reports of air leaking into infusion lines through the BD Q-Syte connector.</p> <p>Risk of air embolus to patients.</p>	<ul style="list-style-type: none"> <li>• Identify, quarantine and do not use devices from affected lot numbers.</li> <li>• Contact Becton Dickinson to arrange for credit or replacement product.</li> </ul>
Action by	
All medical and nursing staff.	
CAS deadlines	Contact
<p>Action underway: 25 November 2009</p> <p>Action complete: 02 December 2009</p>	<p><b>Manufacturer</b>            Customer Services            Becton Dickinson            Tel: 01865 781 517            Fax: 01865 781 501            E-mail: <a href="mailto:val_mummary@europe.bd.com">val_mummary@europe.bd.com</a></p>

[Link to full Medical Device Alert](http://www.mhra.gov.uk/Publications/Safetywarnings/MedicalDeviceAlerts/CON062862)

[www.mhra.gov.uk/Publications/Safetywarnings/MedicalDeviceAlerts/CON062862](http://www.mhra.gov.uk/Publications/Safetywarnings/MedicalDeviceAlerts/CON062862)