

# Medical Device Alert

## Action

Ref: MDA/2009/078 Issued: 08 December 2009 at 11:00

### Device

Intra-aortic balloon pump harness safety system supplied for use with the CS300 balloon pump.

Supplied by Datascope, now Maquet.

Problem	Action
<p>The harness was not designed to be used with the CS300 balloon pump and therefore does not adequately secure it during ambulance transportation.</p>	<ul style="list-style-type: none"> <li>Identify harnesses.</li> <li>Do not use them with the CS300 intra-aortic balloon pump.</li> <li>Return harnesses supplied with the CS300 intra-aortic balloon pump to the manufacturer for a refund and information on alternatives.</li> </ul>
Action by	
<p>All ambulance, medical, nursing and technical staff associated with the transportation of CS300 balloon pumps.</p>	
CAS deadlines	Contact
<p>Action underway: 21 December 2009 Action complete: 04 January 2010</p>	<p><b>Manufacturer</b> Mr Colin Moralee Maquet Ltd Tel: 0191 519 6200 Fax: 0191 519 6201 E-mail: <a href="mailto:cmoralee@maquet.co.uk">cmoralee@maquet.co.uk</a></p>

### Link to full Medical Device Alert

<http://www.mhra.gov.uk/Publications/Safetywarnings/MedicalDeviceAlerts/CON065435>