

# Medical Device Alert

## Action

Ref: MDA/2009/053 Issued: 05 August 2009 at 14:00

### Device

<p>Manual wheelchairs.</p> <p>Breezy 100 and Breezy 300 ranges.</p> <p>Manufactured by Sunrise Medical Ltd.</p>			
	Breezy 110	Breezy 115	Breezy 121
			
	Breezy 141	Breezy 310	Breezy 315

### Problem

These wheelchairs are not suitable for use as a seat in a motor vehicle unless they have been fitted with a transit kit. Without this kit the wheelchair could break loose in a vehicle impact, leading to injury of the wheelchair user or others.

### Action by

All those involved in the provision, prescription, maintenance and use of these wheelchairs. In particular: wheelchair service managers, rehabilitation engineers, maintenance staff and contractors, occupational therapists and patient transport service managers.

### CAS deadlines

Action underway: 07 September 2009  
Action complete: 05 November 2009

### Action

- 1) Trace wheelchairs from the Breezy 100 range manufactured before 24/07/2008 and those from the 300 range manufactured before 28/07/2008.
- 2) Ensure that users and carers are aware that these wheelchairs are not suitable for use as a seat in a vehicle until the transit kit has been fitted.
- 3) Contact Sunrise Medical Ltd to obtain transit kits and fitting information.
- 4) Arrange for the kit to be fitted.
- 5) Arrange for users and transport providers to receive copies of the latest instructions for use.

### Contact

#### Manufacturer

Sunrise Medical Ltd

Customer Services

Tel: 01384 446 666

E-mail: [enquiries@sunmed.co.uk](mailto:enquiries@sunmed.co.uk)

### Link to full Medical Device Alert

<http://www.mhra.gov.uk/Publications/Safetywarnings/MedicalDeviceAlerts/CON054418>